



PAX Übungspuppe Erste-Hilfe PAX First Aid Practice Manikin

Die PAX Übungspuppe Erste-Hilfe wurde von uns in Zusammenarbeit mit aktiven Ausbildern für Erste-Hilfe o.ä. Kurse entwickelt. Diese Puppe eignet sich für das Üben von stabiler Seitenlage und Rettung aus dem Gefahrenbereich. Die Puppe weist ein Eigengewicht von ca. 30 kg auf, was zwar wesentlich leichter als ein Mensch ist, sich aber in Erste-Hilfe-Kursen als ausreichend realitätsnah bewiesen hat. Die Überstreckung des Kopfes und Öffnung des Mundes können ebenfalls geübt werden.

The PAX First Aid Practice manikin was developed by us in cooperation with active trainers for first aid or similar courses. This manikin is suitable for practising stable lateral position and rescue from the danger zone. The manikin has a dead weight of approx. 30 kg, which is considerably lighter than a human being, but has proven to be sufficiently realistic in first aid courses. The hyperextension of the head and opening of the mouth can also be practised.



Von
AUSBILDERN
für die
AUSBILDUNG

From
TRAINERS
for
TRAINING

▲ Einfache Reinigung Easy cleaning

Wischdesinfektion der Oberflächen nach jeder Schulungseinheit
Wipe disinfection of surfaces after each training session

Art.-Nr. / Item No. 272500309

PAX Übungspuppe Erste-Hilfe PAX First Aid Practice Manikin

180cm

30kg

08

TEC



PAX Übungspuppe Erste-Hilfe PAX First Aid Practice Manikin

An den Hacken der Übungspuppe wurde ein austauschbares, aber abriebfestes Material angebracht, um die Übungspuppe vor Beschädigungen durch das Schleifen auf dem Boden zu schützen. Die Oberflächen der Puppe lassen sich einfach mittels Wischdesinfektion nach jedem Kursteilnehmer reinigen. Diese Übungspuppe passt in unseren PAX Bekleidungstrolley XL.

A replaceable but abrasion-resistant material has been attached to the heels of the manikin to protect it from damage caused by dragging on the floor. The surfaces of the manikin can be easily cleaned by wiping disinfection after each course participant. The manikin fits into the PAX Clothing Trolley XL.



PAX Bekleidungstrolley XL / PAX Clothing Trolley XL

