



PAX Strecher Mattress - TMS with temperature management system

The PAX Strecher Mattress - TMS was developed to make lying transport as comfortable as possible for patients. It is attached to the stretcher and held in place by the patient belt system.

The heating system can be switched on long before reaching the destination and offers a pleasant added value, not only during the cold season.




Dimensions :
 Height: approx. 10 cm
 Width: approx. 130 - 196 cm
 Depth: approx. 42 - 55 cm
 (Depending on model.)
Colour: Black
 **Area of application:** Professional

PAX Strecher Mattress - TMS
Sample model images



PAX Strecher Mattress - TMS (Versions)

Item number	Item description
291707203	PAX Strecher Mattress - TMS - Spencer QMX 777+
291717203	PAX Strecher Mattress - TMS - Spencer QMX 779
290777203	PAX Strecher Mattress - TMS - Stollenwerk
290787203	PAX Strecher Mattress - TMS - Stollenwerk CPS
290797203	PAX Strecher Mattress - TMS - Stryker M1
290807203	PAX Strecher Mattress - TMS - Stryker Power Pro
290817203	PAX Strecher Mattress - TMS - Stryker Power Pro XT
290837203	PAX Strecher Mattress - TMS - Ferno Mondial
290847203	PAX Strecher Mattress - TMS - Ferno Viper
290857203	PAX Strecher Mattress - TMS - Kartsana BRAVA SILVER JUPITER
290827203	PAX Strecher Mattress - TMS - Ferno IN-X
290867203	PAX Strecher Mattress - TMS - Kartsana SUPERBRAVA
291443503	PAX Strecher Mattress - TMS - Frestems





it's all about performance

Further information:




PAX TMS Power Pack

The battery for the TMS system

The PAX TMS Power Pack ensures a reliable power supply during use.

Even at -10°C , the battery achieves a runtime of approx. 3 to 3.5 hours, measured with an average 75 kg patient at the highest heating level.
For heavier patients, the runtime is extended accordingly.

The charging time is around 3 hours.
Plug-and-play operation enables simple and safe handling in everyday work.



Item no. 290000000

PAX TMS Power Pack



Item no. 290010000S

PAX TMS Charger 220V - with power cord

PAX TMS charging system

The PAX TMS charging system includes the charging cradle and the power supply (100-240V) with accompanying power cord.

The battery is not included in the scope of delivery.

PAX Vacuum Mattresses

Secure immobilisation
and comfortable patient transport

Also available with
heating function
as TMS variant
for all designs.

