



## PAX RSQ 2-fold ezLock Patient Stretcher

The PAX RSQ 2-fold ezLock Patient Stretcher is a compact and lightweight rescue stretcher designed for professional use in technical rescue operations and in emergency situations where space is limited. The two-fold design combines high stability with a compact build and supports safe patient transport even under demanding conditions.





A key feature is the innovative ezLock locking mechanism:

The one-handed operation reliably prevents the stretcher from folding up accidentally, ensuring a high level of safety and rapid deployment during rescue operations.

Thanks to its compact design, the PAX RSQ 2-fold ezLock patient stretcher can be stowed and transported in the specially designed PAX RSQ Transport Backpack 2-fold ezLock Patient Stretcher.


Item no.: 291234351

### PAX RSQ 2-fold ezLock Patient Stretcher

-  (HxWxD): 183 x 41 x 3 cm
-   (HxWxD): 60,5 x 41 x 10 cm
-  6,00 kg

51

PRO





## PAX RSQ Transport Backpack 2-fold ezLock Patient Stretcher

PAX RSQ Transport Backpack 2-fold ezLock Patient Stretcher was specially developed for the safe transport and space-saving storage of the PAX RSQ 2-fold ezLock Patient Stretcher. Thanks to its foldable design, the transport backpack can be stored in a particularly compact and space-saving manner.

When folded, it can be attached directly to the PAX RSQ 2-fold ezLock Patient Stretcher, ensuring that it is always ready for transport and quickly accessible.

An integrated quick-release mechanism provides rapid access to the patient stretcher and supports efficient action in time-critical emergency situations.

In addition, the backpack is equipped with a winchable handle, making it ideally suited for technical rescue scenarios as well as rope and high-angle rescue operations.



Item no.: 290470151

### PAX RSQ Transport Backpack 2-fold ezLock Patient Stretcher



(HxWxD): 56 x 43 x 12 cm

51



0,73 kg

PRO

