



# Data sheet PAX Roll-Up Stretcher light - tactical

Version according to DIN EN 1865-1

Item No 279576003

## PAX Roll-Up Stretcher light - tactical

Small and lightweight transport device for evacuating people from the danger zone. Due to the reinforcement at the bottom, the PAX Roll-Up Stretcher light - tactical can also be pulled. This allows the product to be operated by only one person. Transport bag with MOLLE. Delivery incl. 1 pc. towing sling with 3 m length.

- 52 cm x Ø14 cm** (rolled up, in bag)
- 212 cm x 50 cm** (unrolled)
- approx. 2 kg** (incl. bag and towing sling)


03

LDPE\*

PRO

\*Material plastic sheet partially recycled.

**Buckle**  
adjustable  
in both directions



**Fixing point**  
for strapping



**Handle**  
2, each at head  
and foot end



**Pouch**  
for towing sling



**lightweight  
and  
compact**



**Rubber loop**  
in lid for optional  
additional strap loop



**Bag with MOLLE**  
loops for universal  
attachment to  
equipment

